

■ Proline® Prospenser

Easy-to-use Bottle-top Dispenser

Proline® Prospenser has been designed for trouble-free and reliable bottle-top dispensing of liquids, including concentrated acids, bases, saline solutions, as well as many organic solvents.

- Dispensing directly from the supply bottle
- Fully autoclavable at 121°C, 2 bar, 20 minutes
- Accuracy of delivery within $\pm 0.5\%$
- Chemically resistant fluid pathway
- Anti-drip closure cap included
- Easy-to-use volume adjustment for reproducible dispensings
- Easy to dismantle for cleaning and maintenance
- Wide range of adapters included to fit the most common bottle sizes
- Optional flexible dispense tube extension (max length 800 mm) with safety handle enables fast and precise dispensing even into narrow tubes
- Each unit is supplied with performance certificate and tested according to ISO 8655

Volume setting easy to adjust

Borosilicate glass barrel protected with a transparent polypropylene sleeve

Anti-drip closure cap

Wide range of bottle adapters included

Chemically resistant liquid pathway



Optional flexible dispense tube extension

- Coiled
- Length 800 mm

