

## □ Ordering Information and Performance Specifications

### Proline® Prospenser

Order Code	Item	Increment	Max Volume	Systematic Error (%)	Random Error (%)
723045	Proline® Prospenser 0.5 – 5 ml (with 25, 28, 32, 38 and 40 mm adaptors)	0,1 ml	5 ml	0,5	0,1
723046	Proline® Prospenser 1 – 10 ml (with 25, 28, 32, 38 and 40 mm adaptors)	0,2 ml	10 ml	0,5	0,1
723047	Proline® Prospenser 2.5 – 25 ml (with 32, 38 and 40 mm adaptors)	0,5 ml	25 ml	0,5	0,1
723048	Proline® Prospenser 5 – 50 ml (with 32, 38 and 40 mm adaptors)	1,0 ml	50 ml	0,5	0,1
721633	Flexible tube extension for 5 and 10 ml Proline® Prospensers				
721634	Flexible tube extension for 25 and 50 ml Proline® Prospensers				

### Prospenser

Order Code	Item	Increment	Max Volume	Systematic Error (%)	Random Error (%)
723049	Prospenser 0.01 – 2.5 ml (with 38, 40 and 45 mm adaptors)	0,05 ml	2,5 ml	0,3	0,1
723050	Prospenser 0.1 – 5 ml (with 38, 40 and 45 mm adaptors)	0,1 ml	5 ml	0,3	0,1
723051	Prospenser 0.2 – 10 ml (with 38, 40 and 45 mm adaptors)	0,2 ml	10 ml	0,3	0,1
723052	Prospenser 1 – 30 ml (with 38, 40 and 45 mm adaptors)	1,0 ml	30 ml	0,3	0,1
723053	Prospenser 1 – 50 ml (with 38, 40 and 45mm adaptors)	1,0 ml	50 ml	0,3	0,1
721998	Dispense tube extension				

### Biotrate Digital Burettes 0 – 30 ml and 0 – 50 ml

Order Code	Item	Increment	Systematic Error (%)	Random Error (%)
723054	Biotrate 0 – 30 ml (with 33, 38 and 45 mm adaptors)	0,01 ml	0,2	0,1
723055	Biotrate 0 – 50 ml (with 33, 38 and 45 mm adaptors)	0,01 ml	0,2	0,1
721998	Expandable delivery jet			