

■ Midi Plus Pipetting Controller

Excellent Performance and Ergonomics

The Midi Plus is a lightweight electronic cordless pipetting controller, which allows aspiration from bottles and tubes, without the arm and hand elevations required in the case of serological or volumetric pipettes. It fits all commonly used 1 - 100 ml glass or plastic pipettes, but can also be used with Sartorius 5 ml and 10 ml disposable tips. The speed can be fine-tuned by applying varying finger pressure to the operating buttons.

The Midi Plus is ideal, for example, in microbiological work: dispensing

into a culture media dish can be performed carefully, drop by drop, without breaking the fine surface of the medium.

- Stepless speed control
- Hydrophobic autoclavable filter protects the device in case of over-aspiration
- Convenient fold-out bench stand supports the unit and pipette when not in use*
- Low battery warning



Stepless speed control

- Allows the aspiration and dispensing speed to be adjusted easily and precisely to suit either large or small volume pipettes.

Fold-out stand*

- Prevents pipette contamination
- Saves space
- Allows the unit to be rested on a table with a pipette attached

Hydrophobic autoclavable filter prevents over-aspiration

Autoclavable pipette adapter set

Fits all commonly used 1-100 ml glass or plastic pipettes

Operating buttons



*Not available in units sold in the United States

□ Features

Pipette types	Plastic or glass 1– 100 ml Pasteur pipettes 5 ml and 10 ml Sartorius pipette tips
Rechargeable during use	Yes
Speed control	Stepless adjustable control
Gravity dispensing	Yes
Stand	Attached support
Weight	207 g
Low battery indicator	Yes
Autoclavable parts	Nose cone, pipette holder and filter

□ Ordering Information

Order Code	Item	Qty/Unit
710931	Midi Plus Pipetting Controller with Universal Charger ¹	1
LH-7129120	Replacement filter, 0,45 µm, non-sterile, autoclavable	5
LH-7129130	Replacement filter, 0,2 µm, sterile	5
LH-711019	Adapter set (standard), autoclavable, including nose-cone and silicone adapter	1
LH-711017	Adapter set for 5 ml pipette tip, autoclavable, including nose-cone and silicone adapter	1
LH-711018	Adapter set for 10 ml pipette tip, autoclavable, including nose-cone and silicone adapter	1
780300	Optifit Tip 5 ml (length 150 mm)	100
780308	Optifit Tip 5 ml (length 150 mm)	1000
780310	Midi Plus Tip 10 ml (length 242 mm)	100

1) Supplied with a universal charger (incl. EU, UK, US | JPN, AUS, KOR and CHN plugs)

